

SYNAPSE

DENTAL PAIN ERASER

HOW TO USE AT HOME



TURN ON

Press button to turn on light.

BRACES PAIN?

Move tips on gums (like painting a pole) for ~20-30 sec per tooth.



TOOTH PAIN?

Move tips ~2 min on gums (like painting a pole - see above) and 10 sec on tooth.

MOUTH SORES?

Full circle AROUND sore for ~30 sec.

NO tips ON sore.



TURN OFF

Press button 2 sec.

Light flashes then off.

Auto shut off in 5 min.

CLEAN

Rinse and dry tips.

Replace cap.



LEARN
MORE



SYNAPSE

DENTAL PAIN ERASER

TIPS FOR USING



USE UNTIL PAIN GOES AWAY
20 sec up to 5 min.

BOTH TIPS MUST TOUCH GUMS OR TOOTH!

Use with a mirror or while being still to make sure (watching TV, reading, etc.)



DON'T SHARE
Use only in your own mouth.

SALIVA HELPS

Too dry or too much saliva? Close lips and swallow.



PREVENT PAIN
Use BEFORE your dental visit, too!

Once you've used the Dental Pain Eraser, it will change your dental experience forever!

